


INFORMATION SHEET



3-PLY HIGH FILTRATION DISPOSABLE MASK

05895808



FEATURES

- Adjustable nose bridge
- Comfortable ear loops
- Easy to adjust

PRODUCT PERFORMANCE

This product is provided in a non-sterile form, consisting of a mask body, nose clip, mask band, rectangular shape.

INSTRUCTIONS FOR USE

Nose bar should be faced outside and the direction of the fold should be downward. Use both hands to expend the mask band and avoid touching the inside of the mask. Put the mask on the face and adjust it to a suitable position.

TECHNICAL DATA

Type	Flat Supra-aural non aseptic type
Size	175×95mm(±5%)
Manufacturer	Zhejiang Yinguang Reflecting Material Manufacturing Co.,Ltd
Tested	Guangzhou Inspection, Testing and Certification Group Co., Ltd.
Efficiency	BFE (Bacterial Filtration Efficiency) > 99% CFU/g (Total Amount Of Bacterial Colony) < 20
Structure	Mask, adjustable nose bridge, ear loops Medial device-type II according to NF EN14683: 2006.
Packaging	50 masks per box

TESTING DATA

BACTERIAL FILTRATION EFFICIENCY	M
COLIFORM GROUP	M
TOTAL AMOUNT OF BACTERIAL COLONY	M
TOTAL AMOUNT OF FUNGUS COLONY	M
PSEUDOMONAS AERUGINOSA	M
STAPHYLOCOCCUS AUREUS	M
HEMOLYTIC STREPTOCOCCUS	M
APPEARANCE[3 PIECES]	M
STRUCTURE AND SIZE[3 PIECES]	M
AIRFLOW RESISTANCE	M

NOTE: "M": MEET THE STANDARD'S REQUIREMENT; "F": FAIL TO MEET THE STANDARD'S REQUIREMENT "---": NO COMMENT